



Flow Cytometry Protocol for Staining of Colonocytes Intracellular Antigen

Protocol:

1. Prepare 25,000 to 50,000 cells (5-8 μ m) per reaction.
2. Wash with 2 ml PBS.
3. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette.
4. Fix cells with 2 ml 1% paraformaldehyde in PBS for 30 minutes.
5. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette.
6. Add 1 ml PBS, suspend them gently by tapping tube.
7. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette.
8. Incubate the cells with 2 ml of 1% BSA, 0.1% Saponin in PBS for 10 minutes
9. Repeat steps 2-3 twice.

For Direct Staining:

10. Add primary antibody with conjugate, usually 5 μ l/50,000 cells (1mg/ml).
11. Incubate the cells at room temperature for 15 minutes or 4 °C for 60 minutes
12. Add 1ml 1% BSA in PBS, suspend them gently by tapping tube.
13. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette.
14. Suspend cells in 1ml 1%BSA in PBS.
15. Analyze cells with flow cytometer.

For Indirect Staining:

10. Add primary antibody without conjugate, usually 5 μ l/50,000 cells (1mg/ml).
11. Incubate the cells at room temperature for 15 minutes or 4 °C for 60 minutes
12. Add 1ml 1% BSA in PBS, suspend them gently by tapping tube.
13. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette. (CONTINUED)



NonInvasive Technologies

8170 Lark Brown Rd. Suite 101, Elkridge, Maryland 21075
www.noninvasivetech.com Tel: 410.799.0551 410.799.1155 Fax: 410.799.0082

14. Suspend cells in 1ml 1%BSA in PBS.
15. Repeat step 13-14 twice.
16. Adding secondary antibody with conjugate, usually 5ul/50,000 cells (1mg/ml).
17. Incubate the cells at room temperature for 15 minutes or 4 °C for 60 minutes.
18. Add 1ml 1% BSA in PBS, suspend them gently by tapping tube.
19. Centrifuge at 2000 rpm for 5 minutes, 4 °C, remove the supernatant using pasteur pipette.
20. Suspend cells in 1ml 1%BSA in PBS.
21. Analyze cells with Flow cytometer.

Flow Cytometer Detectors/Amps setting:

Run FASCComp with Lysis/No wash to reduce signal from debris.

<u>Detectors</u>	<u>Voltage</u>	<u>Mode</u>
FSC	E-0	Log
SSC	280-360	Log
FL1	425-475	Log
FL2	450-500	Log

Notes:

1. First choice of antibody conjugate is PE, ideal for low expression proteins. Second choice is the FITC. Third choice is PerCP, PE-Cy5 or APC.
2. Colonocytes have cell size range from 2 to 8 μm , hence the FSC has 16 levels of differences. It is very difficult to use Quadrant gate to get data. A Regulate gate is a better choice.